



PRESSALIT®

Product Item number	Mounting kit V8653 (4 pcs), for aerated concrete V8653
Product description	EAN:5708590351414 Mounting kit, 4 bolts M8 x 35, plugs FTP M 8, for aerated concrete. Used for: PLUS wall tracks (vertical): R9871, R9872 PLUS sink brackets: R4650, R4750, R4751 PLUS shower seats: R7634, R7635 Mounting kit, 4 bolts M8 x 35, plugs FTP M 8, for aerated concrete
Materials	Stainless steel and plastics

pressalit.com