



PRESSALIT®


Institutional hinge

Item number

D92999
VVS: 616810934
EAN: 5708590300474

Description

Institutional hinge, stainless steel, long bolt (3.9")
Institutional hinge, stainless steel, long bolt (3.9")

**Product capacity****Weight**

0,298 kg

Gross weight

0,298 kg

Weight per crate

0 kg

Qty per crate

0

Qty per pallet

0

This product should be cleaned with hot water or a mild cleaning solution. Do not use chlorine based or corrosive cleaners as it may result in the surface being damaged.

Subject to changes in the product range, technical changes and printing errors.

pressalit.com

PRESSALIT®