




PRESSALIT®

Product Item number	Pressalit Public Classic Waste bin 34 l, brushed steel PY30910
Product description	VVS:777744315 NRF:1438608 RSK:8920394 EAN:5708590392363  Waste bin made of robust stainless steel. Stainless steel body with rounded edges. Equipped with removable, washable polypropylene fabric bag. Capacity approx. 34 liters. Installation can be done with glue or with two screws. There is also the option of mounting with extra anti-theft protection in the back wall. Incl. mounting kit.
Materials	European stainless steel (AISI 304)
Warranty Cleaning	A 2-year warranty is provided for material and manufacturing defects. The product is to be cleaned with hot water or a mild soapy cleaning agent. Do not use chlorine based or corrosive cleaners as it may result in the surface being damaged. Subject to changes in the product range, technical changes and printing errors.

pressalit.com