



United Kingdom Conformity Assessment

United Kingdom RP:
Pressalit Ltd.
100, Longwater Avenue
Green Park
Reading RG2 6GP
United Kingdom
T: 01180923 2696
uk@pressalit.com

This product is produced in compliance with all relevant requirements of:
- Regulation (EU) 2017/745 on medical devices (Risk class I, rule 1)
- Council Directive 2011/65/EC on RoHS

Medical devices - Application of risk management and quality management in design and production

ISO 14971:2019
ISO 13485:2016

The products are manufactured in accordance with the following standards:

ISO 17966:2016

Ry, 28. juli 2026

Mette Elstrøm Kristoffersen, QA/RA Specialist, Pressalit A/S
Responsible

Manufacturer: Pressalit A/S - Pressalitvej 1 DK-8580 Ry - T: +45 8788 8788 - pressalit@pressalit.com

PLUS



R372585108

EAN (GTIN): 5708590368672

Nomenclature (GMDN): 35585

Nomenclature (EMDN): Y181809

Name

EN PLUS support arm with mounting profile, 850 mm, incl. operating panel for TL3/TL4, dark blue. Height adjustable, left hand operated, for fitting to TL3/TL4 (right)

Intended purpose

EN The support arm is used as a suport for a user, when he or she sits down or gets up from a toilet.
Expected service life is 10 years, provided that maintenance is carried out as described in the user manual.