

EAN: 6414150087707  
[static.oras.com/5615](http://static.oras.com/5615)



Oras Clinica washbasin faucet. Firm loop-lever offers good usability to all users and situations. Swivel spout 120° (can be limited to 40°-80°) with nickel-free water channels and laminar water flow. In-built limitation options for water temperature and maximum flow-rate. Connection pipes Cu  $\phi$ 10 mm.

- HealthCare
- Single-lever
- Deck mounted
- Chrome
- Swivel spout, Swivel range limiting option
- Laminar stream
- Single operating lever/handle, Loop shape, Long shape, Hot/Cold symbols
- Limitation option for maximum temperature and flow-rate
- $\phi$  40 mm ceramic cartridge for flow and temperature control
- Copper inlet pipes
- Inner body made of DZR brass, Waterways without nickel coating, 3S-installation system for safe and simple mounting

## Technical data

### Flow properties

Flow-rate at 300 kPa	<b>0.2 l/s</b>
Pressure loss with flow (0.1 l/s)	<b>70 kPa</b>

### Technical properties


DN-size (nominal dimension)	<b>DN15</b>
Hot water supply	<b>max. +80°C</b>
Working pressure	<b>50 - 1000 kPa</b>
Projection	<b>133 mm</b>
Backflow prevention (EN1717)	<b>AA</b>
Material	<b>Brass</b>



### Regulations

EN standard	<b>EN 817</b>
Noise class	<b>I (ISO 3822)</b>


### Approvals and Declarations

STF certificate	<b>EUF129-21002860-TH</b>
KIWA SE	<b>1608</b>
EPD	<b>S-P-06391</b>





### Spare parts



### SP5613

Name	Code
01 Long lever, L=162	1004914V
01 Lever, L=136	1005294V
02 Covering cap	602589V
03 Tightening nut	158726
04 Cartridge	158890
05 Spout, L=132	602520V
06 Spout overturning stopper, 40°-80°	1002207/10
06 Spout overturning stopper, 120°	1001406/10
07 Fixing screw, M8x1, L=130	1000780V
08 Fixing set	602610V
09 Aerator, M24x1 STD PL HC A Laminar	602283V
10 Bearing ring	601097/10
11 X-sealing set	158735
12 Bearing ring	159423/10
13 Headpart, G1/2	602600V
14 Handle for washing machine valve	1005697V
15 Bottom seal	1002371V
16 Supporting plate, 100x60mm	158825