




PRESSALIT®

Product Item number	Pressalit Public Premium Hygienic waste bin 4 l, polished steel PX10925
Product description	VVS: 777743440 NRF: 1438636 RSK: 8920363 EAN: 5708590392127  Wall-mounted sanitary bin with bag dispenser integrated on the side, removable 4 liter sanitary bin and hinged lid. Includes mounting kit and a box with plastic sanitary bags.
Materials	European stainless steel (AISI 304)
Warranty Cleaning	A 2-year warranty is provided for material and manufacturing defects. The product is to be cleaned with hot water or a mild soapy cleaning agent. Do not use chlorine based or corrosive cleaners as it may result in the surface being damaged. Subject to changes in the product range, technical changes and printing errors.

pressalit.com