



PRESSALIT®

Product	Pressalit Trend wall hung toilet with circular flush
Item number	220366
Product description	EAN:5708590398334 The rimless, wall-hung toilet features the efficient Circular Flush system, ensuring a quiet and splash-free flush while keeping everyday cleaning simple.
Product capacity	Max load: 240 kg
Materials	Porcelain
Warranty	Pressalit makes a commitment to remedy material and manufacturing defects for a period of five years from the invoicing date. Subject to changes in the product range, technical changes and printing errors.

pressalit.com