




PRESSALIT®

Product Item number	Pressalit Choice Soapdispenser, 330 ml, brushed steel Q30700
Product description	<p>VVS:774681100 NRF:6190957 RSK:8507854 EAN:5708590356228</p>  <p>Let the Pressalit Choice bathroom accessory series elevate the décor of your bathroom and create coherence with a timeless look. You get a complete range of accessories to suit any bathroom.</p> <p>Elegant soap dispenser for wall mounting. Opens up extra space on the washbasin.</p>
Materials	Stainless steel AISI 304, ABS
Warranty	A 2-year warranty is provided for material and manufacturing defects.
Cleaning	<p>The product is to be cleaned with hot water or a mild soapy cleaning agent. Do not use chlorine based or corrosive cleaners as it may result in the surface being damaged.</p> <p>Subject to changes in the product range, technical changes and printing errors.</p>

pressalit.com