



United Kingdom Conformity Assessment

United Kingdom RP:
 Pressalit Ltd.
 100, Longwater Avenue
 Green Park
 Reading RG2 6GP
 United Kingdom
 T: 01180923 2696
 uk@pressalit.com

This product is produced in compliance with all relevant requirements of:

- Regulation (EU) 2017/745 on medical devices (Risk class I, rule 1)
- Council Directive 2006/42/EC on machinery
- Council Directive 2011/65/EC on RoHS

For the UK:

- Medical Device Regulation (MDR) 2002 on medical devices
- Supply of Machinery (Safety) Regulations 2008
- Restriction of the Use of Certain Hazardous Substances in Electrical and Electronic Equipment Regulations 2012

Medical devices - Application of risk management and quality management in design and production

ISO 14971:2019
 ISO 13485:2016

The products are manufactured in accordance with the following standards:

IEC 60601-1:2005
 ISO 17966:2016

Ry, 23. juni 2026

Mette Elstrøm Kristoffersen, QA/RA Specialist, Pressalit A/S
 Responsible

Manufacturer: Pressalit A/S - Pressalitvej 1 DK-8580 Ry - T: +45 8788 8788 - pressalit@pressalit.com

SELECT TL5



R9211000

EAN (GTIN): 5708590390956

Nomenclature (GMDN): 63055

Nomenclature (EMDN): Y091299

Name

EN SELECT TL5 toilet lifter with tilt function, electrical, for floor outlet

Intended purpose

EN Toilet lifters ensure optimum ergonomics at a toilet, for both seated and standing users with disabilities.
 Expected service life is 10 years, provided that maintenance is carried out as described in the user manual.