



PRESSALIT®

Product Item number	TMA 2 arm support for mounting on toilet, wider fit R1190112000
Colour	VVS:615724200 NRF:6190324 EAN:5708590376233 HMI: 130724 FDA: ILS/D135866 White
Product description	TMA arm supports comes in various widths, and models with or without toilet seat. TMA 2 Arm supports for mounting on toilet, incl. R37 Dania toilet seat. With BL1999 hinge. Dark grey handles.
Adjustability	Fixed mounting
Materials	Powder-coated stainless steel. Soft polyurethane. <i>Flammability-tested (ISO 8191-1). Test report available.</i>
Ergonomy	Toilet seat and cover in thermoset plastic (urea). The folding arms give a good support and the moulded handles are comfortable for the hands and arms.
Features	The support arms can be folded individually.
Warranty	Pressalit makes a commitment to remedy material and manufacturing defects for a period of five years from the invoicing date.
Cleaning	GENERAL The product is to be cleaned with hot water or a mild soapy cleaning agent. Do not use chlorine based, corrosive or abrasive cleaners as it may result in the surface being damaged. The product must be wiped down with a soft rag wrung out in clear water immediately after cleaning, as some materials may otherwise be discolored and stained, and the surface coating may deteriorate.



PRESSALIT®

Product	TMA 2 arm support for mounting on toilet, wider fit
Item number	R1190112000

pressalit.com

Pressalit A/S
Pressalitvej 1
DK-8680 Ry
Tel.: +45 8788 8788
Fax: +45 8788 8789
E-mail: pressalit@pressalit.com

PRESSALIT®