




PRESSALIT®

Product Item number	Diaper basket, 402 x 212 x 174 mm RT767
Product description	VVS:615801767 NRF:7501034 EAN:5708590171074 HMI: 124570 NBS:N13:449  Basket for nappies, A:402 mm x B:212 mm x C: 174 mm. Basket for nappies, A:402 mm x B:212 mm x C: 174 mm.
Materials	Basic material: Steel Surface material: Rilsan
Warranty	Pressalit makes a commitment to remedy material and manufacturing defects for a period of five years from the invoicing date.



PRESSALIT®

Product	Diaper basket, 402 x 212 x 174 mm
Item number	RT767
Cleaning	Baskets The product is to be cleaned with hot water or a mild soapy cleaning agent. Do not use corrosive or abrasive cleaners as it may result in the surface being damaged. For surface disinfection, chlorine-based disinfectants may be used according to the product specification and safety precautions for these products. The following cleaning agents have been tested by Pressalit and may be used, as long as they are diluted according to the instructions on the container. DISINFECTANTS: Clinimax Difficil-S Erisan OXY+ Sprit (husholdningsprit) Suma Bac D10 Virkon Rodalon Desinfektionsspray Lysoformin Spezial CaviCide Sufanios Premium Tvättsprit 70% WipeClean UniversalCleaning Chlorine Disinfection Wipes Klorilli (Kiilto Clean) NON-DISINFECTANTS: Suma Star free D1 Taski 100 Taski 200 Milizid Citro The product must be wiped down with a soft rag wrung out in clear water immediately after cleaning, as some materials may otherwise be discolored and stained, and the surface coating may deteriorate.



PRESSALIT®

Product	Diaper basket, 402 x 212 x 174 mm
Item number	RT767

pressalit.com

Pressalit A/S
Pressalitvej 1
DK-8680 Ry
Tel.: +45 8788 8788
Fax: +45 8788 8789
E-mail: pressalit@pressalit.com

PRESSALIT®