



PRESSALIT®

Product Item number	Pressalit Choice Toilet brush for wall, brushed brass Q70850
Product description	VVS: 778497205 NRF: 7001083 RSK: 8714275 EAN: 5708590372242 DB: 2173273 Let the Pressalit Choice bathroom accessory series elevate the décor of your bathroom and create coherence with a timeless look. You get a complete range of accessories to suit any bathroom. Elegant toilet brush for wall mounting. Opens up extra space on the floor.
Materials	Stainless steel AISI 304, mounting knot in brass
Warranty Cleaning	A 2-year warranty is provided for material and manufacturing defects. This product should be cleaned with hot water or a mild cleaning solution. Do not use chlorine based or corrosive cleaners as it may result in the surface being damaged. Subject to changes in the product range, technical changes and printing errors.

pressalit.com