




PRESSALIT®

Product Item number	Pressalit Public Premium Set with waste bin, paper dispenser and soap dispenser, brushed steel PX30990
Product description	<p> VVS: 773651315 NRF: 1438638 RSK: 8920373 EAN: 5708590392202 </p>  <p> Wall mounting kit consisting of a one-hand operated 200 ml soap dispenser, paper dispenser for approx. 300 paper napkins and a 15 liter waste bin in brushed stainless steel. For wall mounting using glue or screws. The waste bin has an internal holder for standard waste sacks. The soap dispenser has a 200 ml container designed for regular liquid hand soap. The kit has a spring clip lock. Refilling of the paper dispenser is done from the top. </p>
Materials	European stainless steel (AISI 304)
Warranty Cleaning	<p> A 2-year warranty is provided for material and manufacturing defects. </p> <p> The product is to be cleaned with hot water or a mild soapy cleaning agent. Do not use chlorine based or corrosive cleaners as it may result in the surface being damaged. </p> <p> Subject to changes in the product range, technical changes and printing errors. </p>

pressalit.com