



<b>Product</b> <b>Item number</b>	Pressalit Public Premium Toilet paper holder with space for 5 rolls of toilet paper, matte black PX40220
<b>Product description</b>	VVS: 776130531 NRF: 1438653 RSK: 8920375 EAN: 5708590392707  Stainless steel toilet paper holder with space for storing 5 rolls. For wall mounting. Designed for one standard toilet roll and four hidden rolls that automatically drop down. Including mounting kit.  Inspection opening to indicate filling level. Refill through opening in the top.
<b>Materials</b>	European stainless steel (AISI 304)
<b>Warranty</b> <b>Cleaning</b>	A 2-year warranty is provided for material and manufacturing defects.  This product should be cleaned with hot water or a mild cleaning solution. Do not use chlorine based or corrosive cleaners as it may result in the surface being damaged.  Subject to changes in the product range, technical changes and printing errors.