



PRESSALIT®

Institutional hinge

Item number


D92999

VVS: 616810934

EAN: 5708590300474

Description

Institutional hinge, stainless steel, long bolt (100 mm)

**Product capacity****Weight**

0,298 kg

Gross weight

0,298 kg

Weight per crate

0 kg

Qty per crate

0

Qty per pallet

0

The product is to be cleaned with hot water or a mild soapy cleaning agent. Do not use chlorine based or corrosive cleaners as it may result in the surface being damaged.

Subject to changes in the product range, technical changes and printing errors.

pressalit.com

PRESSALIT®