




<b>Product</b> <b>Item number</b>	Pressalit Public Premium Soap dispenser, 200 ml, brushed steel PX30710
<b>Product description</b>	<p>VVS:774684215 NRF:1438645 RSK:8920370 EAN:5708590392172</p>  <p>Stainless steel soap dispenser for wall mounting. One-hand operated and with drip-free soap pump with a 200 ml container. Designed for regular liquid soap or hand lotion. Incl. mounting kit.</p>
<b>Materials</b>	European stainless steel (AISI 316L)
<b>Warranty</b> <b>Cleaning</b>	<p>A 2-year warranty is provided for material and manufacturing defects.</p> <p>The product is to be cleaned with hot water or a mild soapy cleaning agent. Do not use chlorine based or corrosive cleaners as it may result in the surface being damaged.</p> <p>Subject to changes in the product range, technical changes and printing errors.</p>