




PRESSALIT®

Product Item number	Mounting kit V8663 (8 pcs), for aerated concrete V8663
Product description	EAN: 5708590351315  Mounting kit, 8 coach screws Ø10 x 70, plugs GB14, with grey screw covers, for aerated concrete. Used for: SCT 3000 shower change tables: R8435313, R8435318 Mounting kit, 8 coach screws Ø10 x 70, plugs GB14, with grey screw covers, for aerated concrete
Materials	Stainless steel and plastic

pressalit.com

PRESSALIT®