




Product Item number	Pressalit Public Premium Toilet paper holder with space for 5 rolls of toilet paper, polished steel PX10220
Product description	<p>VVS: 776130540 NRF: 1438654 RSK: 8920359 EAN: 5708590392653</p>  <p>Stainless steel toilet paper holder with space for storing 5 rolls. For wall mounting. Designed for one standard toilet roll and four hidden rolls that automatically drop down. Including mounting kit.</p> <p>Inspection opening to indicate filling level. Refill through opening in the top.</p>
Materials	European stainless steel (AISI 304)
Warranty Cleaning	<p>A 2-year warranty is provided for material and manufacturing defects.</p> <p>The product is to be cleaned with hot water or a mild soapy cleaning agent. Do not use chlorine based or corrosive cleaners as it may result in the surface being damaged.</p> <p>Subject to changes in the product range, technical changes and printing errors.</p>