




PRESSALIT®

Product Item number	Pressalit Public Classic Paper dispenser, manual, 500 sheets, matt black PY40110
Product description	VVS: 776135531 NRF: 1438632 RSK: 8920396 EAN: 5708590392387  Stainless steel paper dispenser for wall mounting. For approx. 500 paper napkins. Key-operated opening. Incl. mounting kit. Stainless steel housing. Fills at the top and has an inspection opening to indicate fill level. Recommended paper size (folded): 240-250 mm x 100-115 mm.
Materials	European stainless steel (AISI 304)
Warranty Cleaning	A 2-year warranty is provided for material and manufacturing defects. The product is to be cleaned with hot water or a mild soapy cleaning agent. Do not use chlorine based or corrosive cleaners as it may result in the surface being damaged. Subject to changes in the product range, technical changes and printing errors.

pressalit.com