



PRESSALIT®

Product Item number	Changing table mattress 900 mm R8701432
Color	VVS:777630702 RSK:8565126 EAN:5708590339061 HMI: 125330 Medium gray
Product description	Mattresses for daycare changing tables are available in several colors and in lengths that are adapted to the different changing table lengths. Accessories for changing tables Mattress 29.9" x 35.4" x 30", for 35.4" changing table Phthalate-free
Materials	Cover: Techmadin, phthalate-free Padding: Foam
Cleaning	Upholstery (Cover): Techmadin, phthalate-free Daily cleaning with soapy water and subsequent wiping over with a damp cloth, and wash regularly at 60-90°, without fabric softener. Shrinkage after washing may occur. Do not turn inside out before washing. Can be tumble-dried at max. 60°. Can be disinfected with antiseptic wipes. Alcohol % < 70. Do not use disinfectants intended for hard surfaces. Padding: (foam) The foam contains only DINP softeners which are approved by the EU for the use in articles for children. Do not wash the foam. Do not use the mattress without the cover for reasons of hygiene.



PRESSALIT®

Product	Changing table mattress 900 mm
Item number	R8701432

pressalit.com

Pressalit A/S
Pressalitvej 1
DK-8680 Ry
Tel.: +45 8788 8788
Fax: +45 8788 8789
E-mail: pressalit@pressalit.com

PRESSALIT®

24. september 2026