




# PRESSALIT®

---

## 2 stabilising buffers w/screw

---

<b>Item number</b>	T9807
	EAN:5708590200521
<b>Description</b>	2 stabilising buffers w/screw



---

### Product capacity

---

<b>Weight</b>	0,03 kg
<b>Gross weight</b>	0,03 kg
<b>Weight per crate</b>	0 kg
<b>Qty per crate</b>	0
<b>Qty per pallet</b>	0

The product is to be cleaned with hot water or a mild soapy cleaning agent. Do not use chlorine based or corrosive cleaners as it may result in the surface being damaged.

---

[pressalit.com](http://pressalit.com)

---