



PRESSALIT®

Product Item number	Pressalit Choice Toilet paper holder with plate, brushed steel Q30210
Product description	VVS:776482100 NRF:6190948 RSK:8744048 EAN:5708590336992 Let the Pressalit Choice bathroom accessory series elevate the décor of your bathroom and create coherence with a timeless look. You get a complete range of accessories to suit any bathroom.
Materials	Stainless steel AISI 304, mounting knot in brass
Warranty Cleaning	A 2-year warranty is provided for material and manufacturing defects. This product should be cleaned with hot water or a mild cleaning solution. Do not use chlorine based or corrosive cleaners as it may result in the surface being damaged. Subject to changes in the product range, technical changes and printing errors.

pressalit.com

PRESSALIT®