




PRESSALIT®

Product	Pressalit Trend incl. Pressalit Style D2 1206
Item number	CP1206011CF
	VVS:613326600 NRF:1476341 RSK:7822148 EAN:5708590397672
Colour	White
Product description	 The Pressalit Trend Combi Pack with the Style D2 toilet seat brings contemporary design and personal expression together in one refined solution. The rimless, wall-hung toilet features the efficient Circular Flush system, ensuring a quiet and splash-free flush while keeping everyday cleaning simple. Paired with the Pressalit Style D2, the pack gains an extra layer of exclusivity. The seat's slim, elegant profile matches the toilet's minimalist lines and adds a light, modern touch to the bathroom. With interchangeable design emblems, it's easy to align the seat with the rest of your Pressalit Style accessories or simply switch colours whenever you want a fresh look. Each pack includes two sets of design emblems, and additional options are available if you want to customise further. The Pressalit Style D2 comes in six carefully selected colours, giving you the freedom to craft the atmosphere you want. Behind the aesthetics is the durability Pressalit is known for, supported by a 10-year warranty and features like Soft Close, Lift-Off, and Pressalit SteadyFix® for added stability. The result is a stylish and hygienic combination that suits modern bathrooms where personalisation, comfort, and quality all matter.
Product capacity	Load - ring seat 240 kg



PRESSALIT®

Product	Pressalit Trend incl. Pressalit Style D2 1206
Item number	CP1206011CF
Materials	BOWL: Porcelain TOILET SEAT: The material is colour ingrained duroplast (UF A 10 = urea formaldehyde) that contains no environmentally hazardous substances. UF plastic comprises approx. 70% urea formaldehyde resin, 30% cellulose, as well as 5% minerals, pigments, lubricants and moisture content. Free formaldehyde is released during the production process. Pressalit toilet seats have been measured after curing ($< 0.01 \text{ mg/m}^3$) to emit significantly below the permitted EU limit of 0.08 mg/m^3. ANTIBACTERIAL: This toilet seat features an antibacterial additive with biocidal silver phosphate glass, which helps prevent bacterial growth. BUFFERS: EVA (copolymer made of ethylene and vinyl acetate). DAMPER KIT: Hydraulic damper with damper housing in stainless steel or thermoplastic, mounted with plastic parts in thermoplastic. MOUNT/HINGES: Stainless steel and thermoplastic.
Warranty	Pressalit provides a 10-year warranty on toilet seats and a 5-year warranty on toilet bowls against presentation of a valid receipt. In the event of a complaint, please contact the retailer or store where the product was purchased. The warranty covers defective materials and/or manufacturing defects and/or functional limitations that may occur following proper use and cleaning of the product. Pressalit reserves the right to examine the product at our factory. In the case of a valid complaint, Pressalit will either:
- supply spare parts, or
- replace the product with an identical or a similar model.
If the color of the toilet seat is no longer available, Pressalit reserves the right to supply a similar model in white. Upon replacement or delivery of spare parts, the warranty period will not be extended nor renewed. If replacement is not possible, Pressalit is entitled instead to refund the purchase price. Subject to changes in the product range, technical changes and printing errors.



PRESSALIT®

Product	Pressalit Trend incl. Pressalit Style D2 1206
Item number	CP1206011CF

pressalit.com

Pressalit A/S
Pressalitvej 1
DK-8680 Ry
Tel.: +45 8788 8788
Fax: +45 8788 8789
E-mail: pressalit@pressalit.com

PRESSALIT®

9. september 2026