




PRESSALIT®

Raised buffer with screw 100 mm

Item number	T9804 EAN:5708590183930 Raised buffer with screw 100 mm
Description	Raised buffer with screw 100 mm
Product capacity	
Weight	0,1 kg
Gross weight	0,1 kg
Weight per crate	0 kg
Qty per crate	0
Qty per pallet	0



This product should be cleaned with hot water or a mild cleaning solution. Do not use chlorine based or corrosive cleaners as it may result in the surface being damaged.

pressalit.com

PRESSALIT®