




PRESSALIT®

Product Item number	Freestanding column 37.2, for fixed height PLUS fold down grab bar R9430
Product description	<p>VVS:615851990 NRF:6190667 EAN:5708590320427 NBS:N13:436</p>  <p>Free-standing column, gray, for fixed height PLUS fold down grab bar. Incl. floor-mounting kit.</p> <p>Free-standing column, gray, for fixed height PLUS support arm. Incl. floor-mounting kit.</p>
Materials	Core material: Steel Surface material: Rilsan
Features	Used with freestanding toilets and when it is not possible to wallmount a grab bar. For fixed height fold down grab bar.
Warranty	Pressalit makes a commitment to remedy material and manufacturing defects for a period of five years from the invoicing date.

pressalit.com

Pressalit A/S
Pressalitvej 1
DK-8680 Ry
Tel.: +45 8788 8788
Fax: +45 8788 8789
E-mail: pressalit@pressalit.com

PRESSALIT®