



Free Sales Certificate

The Danish Medicines Agency hereby certifies that the medical devices specified in the attached list are manufactured by:

Pressalit A/S
Pressalitvej 1
DK-8680 Ry

Medical devices which are CE marked in conformity with Directive 93/42/EEC meet the essential requirements for safety and performance. They may therefore be manufactured and marketed in Denmark and exported without any approval from the Danish Medicines Agency.



Ugur Erman

Ugur Erman

Valid from: **15 March 2018**
Valid Until: **15 March 2020**



Article no. Description

R4950000	SELECT wash basin lifter, floor standing, with white cover plate. Electrically height adjustable with lever control
R4950112	SELECT wash basin lifter, floor standing, with anthracite cover plate. Electrically height adjustable with lever control
R4950277	SELECT wash basin lifter, floor standing, with lime green cover plate. Electrically height adjustable with lever control
R4951000	SELECT wash basin lifter, floor standing, with white cover plate. Electrically height adjustable with wired hand control
R4951112	SELECT wash basin lifter, floor standing, with anthracite cover plate. Electrically height adjustable with wired hand control
R4951277	SELECT wash basin lifter, floor standing, with lime green cover plate. Electrically height adjustable with wired hand control
R4952000	SELECT wash basin lifter, floor standing, with white cover plate. Electrically height adjustable with lever control. Incl. safety stop function
R4952112	SELECT wash basin lifter, floor standing, with anthracite cover plate. Electrically height adjustable with lever control. Incl. safety stop function
R4953000	SELECT wash basin lifter, floor standing, with white cover plate. Electrically height adjustable with wired hand control. Incl. safety stop function
R4953112	SELECT wash basin lifter, floor standing, with anthracite cover plate. Electrically height adjustable with wired hand control. Incl. safety stop function
R8012000	SELECT toilet lifter, manual, floor outlet, 300 mm height adjustable with crank handle, white
R8012112	SELECT toilet lifter, manual, floor outlet, 300 mm height adjustable with crank handle, anthracite
R8012277	SELECT toilet lifter, manual, floor outlet, 300 mm height adjustable with crank handle, lime green
R8022000	SELECT toilet lifter, manual, wall outlet, 300 mm height adjustable with crank handle, white
R8022112	SELECT toilet lifter, manual, wall outlet, 300 mm height adjustable with crank handle, anthracite





R8022277	SELECT toilet lifter, manual, wall outlet, 300 mm height adjustable with crank handle, lime green
R8112000	SELECT toilet lifter, manual, with integrated vertical side profiles, floor outlet, 300 mm height adjustable with crank handle, white
R8112112	SELECT toilet lifter, manual, with integrated vertical side profiles, floor outlet, 300 mm height adjustable with crank handle, anthracite
R8112277	SELECT toilet lifter, manual, with integrated vertical side profiles, floor outlet, 300 mm height adjustable with crank handle, lime green
R8122000	SELECT toilet lifter, manual, with integrated vertical side profiles, wall outlet, 300 mm height adjustable with crank handle, white
R8122112	SELECT toilet lifter, manual, with integrated vertical side profiles, wall outlet, 300 mm height adjustable with crank handle, anthracite
R8122277	SELECT toilet lifter, manual, with integrated vertical side profiles, wall outlet, 300 mm height adjustable with crank handle, lime green
R8212000	SELECT toilet lifter, powered with wired hand control, floor outlet, 400 mm height adjustable with wired remote control, white
R8212112	SELECT toilet lifter, powered with wired hand control, floor outlet, 400 mm height adjustable with wired remote control, anthracite
R8212277	SELECT toilet lifter, powered with wired hand control, floor outlet, 400 mm height adjustable with wired remote control, lime green
R8222000	SELECT toilet lifter, powered with wired hand control, wall outlet, 400 mm height adjustable with wired remote control, white
R8222112	SELECT toilet lifter, powered with wired hand control, wall outlet, 400 mm height adjustable with wired remote control, anthracite
R8222277	SELECT toilet lifter, powered with wired hand control, wall outlet, 400 mm height adjustable with wired remote control, lime green
R8312000	SELECT toilet lifter, powered with wired hand control, with integrated vertical side profiles, floor outlet, 400 mm height adjustable with wired remote control, white
R8312112	SELECT toilet lifter, powered with wired hand control, with integrated vertical side profiles, floor outlet, 400 mm height adjustable with wired remote control, anthracite





R8312277	SELECT toilet lifter, powered with wired hand control, with integrated vertical side profiles, floor outlet, 400 mm height adjustable with wired remote control, lime green
R8322000	SELECT toilet lifter, powered with wired hand control, with integrated vertical side profiles, wall outlet, 400 mm height adjustable with wired remote control, white
R8322112	SELECT toilet lifter, powered with wired hand control, with integrated vertical side profiles, wall outlet, 400 mm height adjustable with wired remote control, anthracite
R8322277	SELECT toilet lifter, powered with wired hand control, with integrated vertical side profiles, wall outlet, 400 mm height adjustable with wired remote control, lime green



LÆGEMIDDELSTYRELSEN
DANISH MEDICINES AGENCY