

IAPMO RESEARCH AND TESTING, INC.

5001 E. Philadelphia Street, Ontario, CA 91761 • Phone (909) 472-4100 • Fax (909) 472-4244 • www.iapmort.org



CERTIFICATE OF LISTING



IAPMO Research and Testing, Inc. is a product certification body which tests and inspects samples taken from the supplier's stock or from the market or a combination of both to verify compliance to the requirements of applicable codes and standards. This activity is coupled with periodic surveillance of the supplier's factory and warehouses as well as the assessment of the supplier's Quality Assurance System. This listing is subject to the conditions set forth in the characteristics below and is not to be construed as any recommendation, assurance or guarantee by IAPMO Research and Testing, Inc. of the product acceptance by Authorities Having Jurisdiction.

Issued To:

PRESSALIT A/S

PRESSALITVEJ 1 RY DK-8680 , DENMARK

Product:

Wall Mounted, Pedestal Mounted, Adjustable, Elevating, Tilting, and Pivoting Lavatory, Sink and Shampoo Bowl Carrier Systems and Drain Waste Systems

Products are in compliance with the following code(s):

N/A - CLASSIFIED

(This is not a listing for code compliance)

Products are in compliance with the following standard(s)

ASME A112.19.12-2014

File Number: C-11505

Revised Date: September 24, 2021

Effective Date: January 2019

Void After: January 2024*


Chairman, Product Certification Committee





Chief Technical Service Officer

*This certificate is not evidence of current listing. To verify listing status, visit the IAPMO R&T Product Listing Directory at pld.iapmo.org

This listing period is based upon the last date of the month indicated on the Effective Date and Void After Date shown above. Any change in material, manufacturing process, marking or design without having first obtained the approval of the Product Certification Committee, or any evidence of non-compliance with applicable codes and standards or of inferior workmanship, may be deemed sufficient cause for revocation of this listing. Production of or reference to this form for advertising purposes may be made only by specific written permission of IAPMO Research and Testing, Inc. Any alteration of this certificate could be grounds for revocation of the listing. This document shall be reproduced in its entirety.

IAPMO RESEARCH AND TESTING, INC.

CERTIFICATE OF LISTING



Issued To: PRESSALIT A/S

Effective Date: January 2019

File Number: C-11505

Revised Date: September 24, 2021

Product: Wall Mounted, Pedestal Mounted, Adjustable, Elevating, Tilting, and Pivoting Lavatory, Sink and Shampoo Bowl Carrier Systems and Drain Waste Systems

Void After: January 2024

This certificate is not evidence of current listing. To verify listing status, visit the IAPMO R&T Product Listing Directory at pld.iapmo.org

Identification:

System shall be permanently marked with the Manufacturer's name or trademark and model number. System shall also bear the Classified certification mark.


Characteristics:

A flexible feed and waste systems to enable the adjustment of the Pressalit Care wash basin brackets. The adjustment can be vertically or vertically as well as horizontally adjusted in the height. To be installed in accordance with the manufacturer's instructions. The flexible feed system has the dimension of 10 mm pipe in one end and 3/8" swivel nut in the other end.

Products listed on this certificate have been tested by an IAPMO R&T recognized laboratory. This recognition has been granted based upon the laboratory's compliance to the applicable requirements of ISO/IEC 17025.

IAPMO RESEARCH AND TESTING, INC.

CERTIFICATE OF LISTING



Issued To: PRESSALIT A/S

File Number: C-11505

Product: Wall Mounted, Pedestal Mounted, Adjustable, Elevating, Tilting, and Pivoting Lavatory, Sink and Shampoo Bowl Carrier Systems and Drain Waste Systems

Effective Date: January 2019

Revised Date: September 24, 2021

Void After: January 2024

This certificate is not evidence of current listing. To verify listing status, visit the IAPMO R&T Product Listing Directory at pld.iapmo.org

Models

Waste Pipe with Flexible Feed Pipes:

Model Number	Description
LA8933	Waste pipe 520mm/flexible feed pipe 700mm
R2084	waste pipe(R2066), with 1200mm Flexible Feed Pipes
R2085	waste pipe(R2063), with 700mm Flexible Feed Pipes
R2087	Waste pipe 700mm/flexible feed pipe 700mm

Brackets/Carriers Note: to be used as a system with listed models

Model Number	Description
R2010	MATRIX supports, pale grey, fixed height, wall-mounted. For MATRIX MEDIUM and ANGLE vanity sinks
R4826	MATRIX MEDIUM sink bracket, wall-mounted, with white installation cover, incl. feed/waste system LA8933. Manually height adjustable with crank handle
R4808	MATRIX MEDIUM sink bracket, wall-mounted, with white installation cover, incl. feed/waste system LA8933. Fixed height
R4828	MATRIX MEDIUM sink bracket, wall-mounted, with white installation cover, incl. feed/waste system LA8933. Electrically height adjustable with operating buttons (in front edge of sink)
R2015	MATRIX supports, pale grey, manually height adjustable, wall-mounted. For MATRIX MEDIUM and ANGLE vanity sinks