




PRESSALIT®

Product Item number	Pressalit Public Premium Toilet paper holder for 1 roll of toilet paper + 1 spare roll, brushed steel PX30210
Product description	VVS: 776130215 NRF: 1438646 RSK: 8920365 EAN: 5708590392660  Toilet paper holder in steel housing with space for spare roll. For wall mounting. The holder is lockable and includes mounting kit.
Materials	European stainless steel (AISI 304)
Warranty Cleaning	A 2-year warranty is provided for material and manufacturing defects. The product is to be cleaned with hot water or a mild soapy cleaning agent. Do not use chlorine based or corrosive cleaners as it may result in the surface being damaged. Subject to changes in the product range, technical changes and printing errors.

pressalit.com