




PRESSALIT®

| | |
|--------------------------------------|--|
| Product Item number | Pressalit Public Classic Paper dispenser, manual, 500 sheets, polished steel PY10110 |
| Product description | <p>VVS:776135540 NRF:1438633 RSK:8920380 EAN:5708590392240</p>  <p>Stainless steel paper dispenser for wall mounting. For approx. 500 paper napkins. Key-operated opening. Incl. mounting kit.</p> <p>Stainless steel housing. Fills at the top and has an inspection opening to indicate fill level.</p> |
| Materials | European stainless steel (AISI 304) |
| Warranty Cleaning | <p>A 2-year warranty is provided for material and manufacturing defects.</p> <p>The product is to be cleaned with hot water or a mild soapy cleaning agent. Do not use chlorine based or corrosive cleaners as it may result in the surface being damaged.</p> <p>Subject to changes in the product range, technical changes and printing errors.</p> |

pressalit.com

PRESSALIT®