




PRESSALIT®

Product Item number	Pressalit Public Premium Toilet paper holder for 1 large roll of toilet paper, brushed steel PX30225
Product description	VVS: 776130115 NRF: 1438649 RSK: 8920367 EAN: 5708590392677  Toilet paper holder in large steel housing with space for rolls of max. Ø 280 x 100 mm. For wall mounting. The holder is lockable and includes mounting kit.
Materials	European stainless steel (AISI 304)
Warranty Cleaning	A 2-year warranty is provided for material and manufacturing defects. The product is to be cleaned with hot water or a mild soapy cleaning agent. Do not use chlorine based or corrosive cleaners as it may result in the surface being damaged. Subject to changes in the product range, technical changes and printing errors.

pressalit.com