



PRESSALIT®

Product Item number	TMA 1 arm support for mounting on toilet R1170000
Color	VVS:615723200 NRF:6191649 RSK:7828180 EAN:5708590211701 HMI: 35270 NBS:N13:436 FDA:ILS/D135866 White
Product description	TMA arm supports comes in various widths, and models with or without toilet seat. TMA 1 arm supports for mounting on toilet, incl. R37 Dania toilet seat. With BL1999 hinge. Dark gray handles.
Adjustability	Fixed mounting
Materials	Powder-coated stainless steel. Soft polyurethane. <i>Flammability-tested (ISO 8191-1). Test report available.</i> Toilet seat and cover in thermoset plastic (urea).
Ergonomy	The folding arms provide firm support and the molded handles are comfortable to hold on to.
Features	The arm supports can be folded individually.
Warranty	Pressalit makes a commitment to remedy material and manufacturing defects for a period of five years from the invoicing date.
Cleaning	This product should be cleaned with hot water or a mild cleaning solution. Do not use chlorine based, corrosive or abrasive cleaners as it may result in the surface being damaged. The product must be wiped down with a soft rag wrung out in clear water immediately after cleaning, as some materials may otherwise be discolored and stained, and the surface coating may deteriorate.



PRESSALIT®

Product	TMA 1 arm support for mounting on toilet
Item number	R1170000

pressalit.com

Pressalit A/S
Pressalitvej 1
DK-8680 Ry
Tel.: +45 8788 8788
Fax: +45 8788 8789
E-mail: pressalit@pressalit.com

PRESSALIT®

24. september 2026