



# PRESSALIT®

|                                      |                                                                                                                                                                                                                                                                                                                                                                                                            |
|--------------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product</b><br><b>Item number</b> | Wall-mounted toilet 27 4/8"<br>R2099000                                                                                                                                                                                                                                                                                                                                                                    |
| <b>Color</b>                         | VVS:612928000<br>NRF:6191026<br>RSK:7822154<br>EAN:5708590345284<br>HMI: 105700<br>White                                                                                                                                                                                                                                                                                                                   |
| <b>Product description</b>           | Toilet, wall-mounted, distance from wall to the front: 27 4/8"<br>Excl. mounting kit.<br>For mounting directly on wall, a mounting kit from V&B must be ordered from Pressalit at an additional cost. For mounting on built-in cistern or Pressalit toilet lifter, mounting kit is included in the delivery.<br><br>Toilet, wall-mounted, distance from wall to the front: 27 4/8".<br>Excl. mounting kit. |
| <b>Product capacity</b>              | Max load: 880 lbs.                                                                                                                                                                                                                                                                                                                                                                                         |
| <b>Materials</b>                     | Sanitary porcelain                                                                                                                                                                                                                                                                                                                                                                                         |
| <b>Ergonomy</b>                      | 27 4/8" is recommended for wheelchair users as it provides sufficient space for a caregiver.<br>Wall mounted toilet makes mounting possible at individual height. Leaves space for foot support, the legs of a floor lift, and the caregiver's feet.                                                                                                                                                       |
| <b>Warranty</b>                      | Pressalit makes a commitment to remedy material and manufacturing defects for a period of five years from the invoicing date.                                                                                                                                                                                                                                                                              |



# PRESSALIT®

|                    |                             |
|--------------------|-----------------------------|
| <b>Product</b>     | Wall-mounted toilet 27 4/8" |
| <b>Item number</b> | R2099000                    |

---

[pressalit.com](http://pressalit.com)

---

Pressalit A/S  
Pressalitvej 1  
DK-8680 Ry  
Tel.: +45 8788 8788  
Fax: +45 8788 8789  
E-mail: [pressalit@pressalit.com](mailto:pressalit@pressalit.com)

# PRESSALIT®