




PRESSALIT®

Product Item number	Pressalit Public Classic Paper dispenser, manual, 500 sheets, matt black PY40110
Product description	<p>VVS:776135531 NRF:1438632 RSK:8920396 EAN:5708590392387</p>  <p>Stainless steel paper dispenser for wall mounting. For approx. 500 paper napkins. Key-operated opening. Incl. mounting kit.</p> <p>Stainless steel housing. Fills at the top and has an inspection opening to indicate fill level.</p> <p>Recommended paper size (folded): 240-250 mm x 100-115 mm.</p>
Materials	European stainless steel (AISI 304)
Warranty Cleaning	<p>A 2-year warranty is provided for material and manufacturing defects.</p> <p>The product is to be cleaned with hot water or a mild soapy cleaning agent. Do not use chlorine based or corrosive cleaners as it may result in the surface being damaged.</p> <p>Subject to changes in the product range, technical changes and printing errors.</p>

pressalit.com