




PRESSALIT®

Product	Mounting kit V8654 (4 pcs), for aerated concrete
Item number	V8654
Product description	EAN:5708590351407
	
Product description	Mounting kit, 4 countersunk bolts M8 x 35, plugs FTP M 8, for aerated concrete.
	Used for: PLUS support arms: R360070, R360085 PLUS folding seat: R5510, R5504 PLUS shower seats: R7304, R7404
Materials	Mounting kit, 4 countersunk bolts M8 x 35, plugs FTP M 8, for aerated concrete
	Stainless steel and plastic

pressalit.com

PRESSALIT®