




<b>Product</b> <b>Item number</b>	Pressalit Public Premium Toilet paper holder with space for 5 rolls of toilet paper, brushed steel PX30220
<b>Product description</b>	<p>VVS: 776130515 NRF: 1438652 RSK: 8920366 EAN: 5708590392141</p>  <p>Stainless steel toilet paper holder with space for storing 5 rolls. For wall mounting. Designed for one standard toilet roll and four hidden rolls that automatically drop down. Including mounting kit.</p> <p>Inspection opening to indicate filling level. Refill through opening in the top.</p>
<b>Materials</b>	European stainless steel (AISI 304)
<b>Warranty</b> <b>Cleaning</b>	<p>A 2-year warranty is provided for material and manufacturing defects.</p> <p>The product is to be cleaned with hot water or a mild soapy cleaning agent. Do not use chlorine based or corrosive cleaners as it may result in the surface being damaged.</p> <p>Subject to changes in the product range, technical changes and printing errors.</p>