




# PRESSALIT®

<b>Product</b> <b>Item number</b>	Pressalit Public Classic Soap dispenser with sensor, 850 ml, brushed steel PY30715
<b>Product description</b>	<p>VVS: 774565215 NRF: 1438639 RSK: 8920389 EAN: 5708590392738</p>  <p>Soap dispenser with stainless steel sensor for wall mounting. Drip-free soap pump with 850 ml container. Designed for regular liquid soap or hand cream. Lockable solution. Incl. mounting kit.</p>
<b>Materials</b>	European stainless steel (AISI 316L)
<b>Warranty</b> <b>Cleaning</b>	<p>A 2-year warranty is provided for material and manufacturing defects.</p> <p>The product is to be cleaned with hot water or a mild soapy cleaning agent. Do not use chlorine based or corrosive cleaners as it may result in the surface being damaged.</p> <p>Subject to changes in the product range, technical changes and printing errors.</p>

[pressalit.com](http://pressalit.com)