



PRESSALIT®

Product Item number	Pressalit Public Classic Paper dispenser, manual, 500 sheets, brushed steel PY30110
Product description	VVS:776135515 NRF:1438631 RSK:8920387 EAN:5708590392318 Stainless steel paper dispenser for wall mounting. For approx. 500 paper napkins. Key-operated opening. Incl. mounting kit. Stainless steel housing. Fills at the top and has an inspection opening to indicate fill level. Recommended paper size (folded): 240-250 mm x 100-115 mm.
Materials	European stainless steel (AISI 304)
Warranty Cleaning	A 2-year warranty is provided for material and manufacturing defects. This product should be cleaned with hot water or a mild cleaning solution. Do not use chlorine based or corrosive cleaners as it may result in the surface being damaged. Subject to changes in the product range, technical changes and printing errors.

pressalit.com