



United Kingdom RP: Pressalit Ltd. 100, Longwater Avenue Green Park Reading RG2 6GP United Kingdom T: 01180923 2696 uk@pressalit.com

The undersigned hereby declares that the below-mentioned products comply with the relevant requirements of the following regulation/directives:

- Medical Device Regulation 2002

- The Supply of Machinery (Safety) Regulations 2008.

	DS/ISO 17966:2016
The products are manufactured in accordance with the following standards:	

Model	R8112000	GMDN	63053	EAN	5708590345352
Description					
SELECT TL2 toilet bracket, for floor outlet, with wide side profiles, white cover plate. Manually height adjustable with					
crank handle					

Model	R8112112	GMDN	63053	EAN	5708590345369	
Description						
SELECT TL2 toilet bracket, for floor outlet, with wide side profiles, anthracite cover plate. Manually height adjustable with						
crank handle						

Model	R8122000	GMDN	63053	EAN	5708590345383	
Description						
SELECT TL2 toilet bracket, for wall outlet, with wide side profiles, white cover plate. Manually height adjustable with						
crank handle						

Model	R8122112	GMDN	63053	EAN	5708590345390
Description					
SELECT TL2 toilet bracket, for wall outlet, with wide side profiles, anthracite cover plate. Manually height adjustable with					
crank handle					

Ry, 5/17/2022

Trine Danielsen, Regulatory Consultant, Pressalit A/S

Signed by the person, who is authorized to compile the technical file, and empowered by the manufacturer to draw up this declaration.



Manufacturer:

Pressalit A/S - Pressalitvej 1 – DK-8680 Ry - T: +45 8788 8788 - pressalit@pressalit.com