




PRESSALIT®

Product Item number	Pressalit Choice Soapdispenser, 330 ml, polished steel Q10700
Product description	VVS:774681140 NRF:6190941 RSK:8507853 EAN:5708590356211  Let the Pressalit Choice bathroom accessory series elevate the décor of your bathroom and create coherence with a timeless look. You get a complete range of accessories to suit any bathroom. Elegant soap dispenser for wall mounting. Opens up extra space on the washbasin.
Materials	Stainless steel AISI 304, ABS
Warranty Cleaning	A 2-year warranty is provided for material and manufacturing defects. The product is to be cleaned with hot water or a mild soapy cleaning agent. Do not use chlorine based or corrosive cleaners as it may result in the surface being damaged. Subject to changes in the product range, technical changes and printing errors.

pressalit.com