




PRESSALIT®

Product	Wall-hung toilet. rimless, 540 mm
Item number	R2098000
	VVS:613232000 NRF:6191063 EAN:5708590377155
Colour	White
Product description	 Toilet, wall-hung, distance from wall to the front: 540 mm Excl. mounting kit. Toilet, wall-hung, distance from wall to the front: 540 mm. Excl. mounting kit.
Product capacity	Max load: 400 kg.
Materials	Ceramic
Ergonomy	Wall-hung toilet makes mounting possible at individual height. Leaves space for foot support, the legs of a mobile hoist, and the carer's feet.
Warranty	Pressalit makes a commitment to remedy material and manufacturing defects for a period of five years from the invoicing date.

pressalit.com

PRESSALIT®