



# PRESSALIT®


## Raised hinge 50 mm

**Item number**

T9802  
EAN:5708590183923

**Description**

Raised hinge 50 mm

**Product capacity****Weight**

0,37 kg

**Gross weight**

0,37 kg

**Weight per crate**

0 kg

**Qty per crate**

0

**Qty per pallet**

0

This product should be cleaned with hot water or a mild cleaning solution. Do not use chlorine based or corrosive cleaners as it may result in the surface being damaged.

[pressalit.com](http://pressalit.com)