



PRESSALIT®

Product Item number	Pressalit Choice Shelf, matt black Q40410
Product description	VVS:772471331 NRF:7001575 RSK:8890000 EAN:5708590381831 Let the Pressalit Choice bathroom accessory series elevate the décor of your bathroom and create coherence with a timeless look. You get a complete range of accessories to suit any bathroom.
Materials	Stainless steel AISI 304
Warranty Cleaning	A 2-year warranty is provided for material and manufacturing defects. This product should be cleaned with hot water or a mild cleaning solution. Do not use chlorine based or corrosive cleaners as it may result in the surface being damaged. Subject to changes in the product range, technical changes and printing errors.

pressalit.com

PRESSALIT®