




<b>Product</b> <b>Item number</b>	Pressalit Choice Towel hook, pin, 2 pcs., matt black Q40150
<b>Product description</b>	<p>VVS:775295231 NRF:6190343 RSK:8754892 EAN:5708590366876</p>  <p>Let the Pressalit Choice bathroom accessory series elevate the décor of your bathroom and create coherence with a timeless look. You get a complete range of accessories to suit any bathroom.</p> <p>Use the towel hook in the bathroom and in the kitchen. That way, you can let the style be consistent throughout the house.</p>
<b>Materials</b>	Stainless steel AISI 304, mounting knot in brass
<b>Warranty</b> <b>Cleaning</b>	<p>A 2-year warranty is provided for material and manufacturing defects.</p> <p>The product is to be cleaned with hot water or a mild soapy cleaning agent. Do not use chlorine based or corrosive cleaners as it may result in the surface being damaged.</p> <p>Subject to changes in the product range, technical changes and printing errors.</p>