



PRESSALIT®

Product Item number	Pressalit Public Classic Disinfection dispenser with sensor, 700 ml, brushed steel PY30725
Product description	VVS: 774565615 NRF: 1438614 RSK: 8920390 EAN: 5708590392745 Disinfection dispenser with stainless steel sensor for wall mounting. One-handed disinfection pump with 700 ml tank. Designed for standard hand sanitizers. Lockable solution. Including mounting kit. Inspection opening indicates filling level.
Materials	European stainless steel (AISI 316L)
Warranty Cleaning WEEE	A 2-year warranty is provided for material and manufacturing defects. This product should be cleaned with hot water or a mild cleaning solution. Do not use chlorine based or corrosive cleaners as it may result in the surface being damaged. Please enquire of your local Pressalit subsidiary or distributor (in countries where there is no subsidiary) that will arrange for recycling of electronic waste in accordance with the WEEE directive. Subject to changes in the product range, technical changes and printing errors.

pressalit.com