




# PRESSALIT®

<b>Product</b> <b>Item number</b>	Pressalit Public Classic Soap dispenser, 1100 ml, brushed steel PY30710
<b>Product description</b>	<p>VVS:774565115 NRF:1438642 RSK:8920388 EAN:5708590392325</p>  <p>Stainless steel soap dispenser for wall mounting. One-handed, drip-free soap pump with 1100 ml container. Designed for regular liquid soap or hand cream. Lockable solution. Incl. mounting kit.</p>
<b>Materials</b>	European stainless steel (AISI 316L)
<b>Warranty</b> <b>Cleaning</b>	<p>A 2-year warranty is provided for material and manufacturing defects.</p> <p>The product is to be cleaned with hot water or a mild soapy cleaning agent. Do not use chlorine based or corrosive cleaners as it may result in the surface being damaged.</p> <p>Subject to changes in the product range, technical changes and printing errors.</p>

[pressalit.com](http://pressalit.com)