




PRESSALIT®

Product	Soap holder
Item number	RT190000
	VVS: 777564500 NRF: 6311237 EAN: 5708590333755 HMI: 88858
Color	White
Product description	Soap holder for grab bars Soap holder for grab bars
Materials	Nylon (PA6)
Warranty	Pressalit makes a commitment to remedy material and manufacturing defects for a period of five years from the invoicing date.



pressalit.com

Pressalit A/S
Pressalitvej 1
DK-8680 Ry
Tel.: +45 8788 8788
Fax: +45 8788 8789
E-mail: pressalit@pressalit.com

PRESSALIT®