



# PRESSALIT®

---

## Splashguard, transparent mat, for 100 mm raised toilet seats

---

|                  |                                                                                                                                  |                                                                                     |
|------------------|----------------------------------------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|
| Item number      | T9515<br>VVS:616810515<br>NRF:6190362<br>EAN:5708590070551                                                                       |  |
| Description      | Splashguard, transparent mat, for 100 mm raised toilet seats<br><br>Splashguard, transparent mat, for 100 mm raised toilet seats |                                                                                     |
|                  |                                                                                                                                  |                                                                                     |
|                  |                                                                                                                                  |                                                                                     |
| Materials        | Polypropylene (PP)                                                                                                               |                                                                                     |
| Product capacity |                                                                                                                                  |                                                                                     |
| Weight           | 0,5 kg                                                                                                                           |                                                                                     |
| Gross weight     | 0,5 kg                                                                                                                           |                                                                                     |
| Weight per crate | 0 kg                                                                                                                             |                                                                                     |
| Qty per crate    | 0                                                                                                                                |                                                                                     |
| Qty per pallet   | 0                                                                                                                                |                                                                                     |

---

[pressalit.com](http://pressalit.com)

---