



PRESSALIT®

Pressalit Style D300 1204

Toilet seat with soft close and lift-off incl. hinge in stainless steel
Incl. design emblems, brushed brass and brushed steel

Features

Lift-off, Soft close

Item number

1204231-DN2999
VVS:614963231
NRF:1587554
RSK:7824153
EAN:5708590399058

Pressalit Style D300 is created for those who want more than just a toilet seat. Its slim, refined profile creates a light and exclusive look that complements and fits any modern bathroom. With interchangeable side design emblems, you can create your bathroom style and a cohesive look with the rest of your Pressalit Style bathroom accessories and tailor the details to match your personal taste and style. The seat includes two sets of design emblems, making it easy to change the colour and appearance whenever you want a refresh. Additional design emblems in other colours and finishes are also available for purchase. Pressalit Style D300 is available in six beautiful colours: white, matt white, black, matt black, matt dusty green and matt taupe. These shades allow you to create exactly the atmosphere you want in your bathroom. Behind the design lies Pressalit's renowned quality, a 10-year warranty, and features such as Soft Close for a quiet, luxurious closing and Lift-off for easy and efficient cleaning. In addition, Pressalit SteadyFix® ensures exceptional stability.






Colour

Matt black

Description

PRESSALIT toilet seat with cover, model "Pressalit Style D300", art. no. 1204 made of colour ingrained duroplast, with built-in soft close function, incl. DN2 universal hinge with lift-off in stainless steel.
- centre distance: 144-224 mm
- load - ring seat 240 kg
- 10 year guarantee





PRESSALIT®

Pressalit Style D300 1204

Toilet seat with soft close and lift-off incl. hinge in stainless steel
Incl. design emblems, brushed brass and brushed steel

Design type	Sandwich
Designer	Pressalit
Design year	2026
Materials	TOILET SEAT: The material is colour ingrained duroplast (UF A 10 = urea formaldehyde) that contains no environmentally hazardous substances. UF plastic comprises approx. 70% urea formaldehyde resin, 30% cellulose, as well as 5% minerals, pigments, lubricants and moisture content. Free formaldehyde is released during the production process. Pressalit toilet seats have been measured after curing ($< 0.01 \text{ mg/m}^3$) to emit significantly below the permitted EU limit of 0.08 mg/m^3. BUFFERS: EVA (copolymer made of ethylene and vinyl acetate). DAMPER KIT: Hydraulic damper with damper housing in stainless steel or thermoplastic, mounted with plastic parts in thermoplastic. MOUNT/HINGES: Stainless steel and thermoplastic.
Product capacity	Load - ring seat 240 kg
Weight	2,748 kg
Carton size	457x389x68 mm
Gross weight	3,172 kg
Crate size	469x401x348 mm
Weight per crate	15,86 kg
Qty per crate	5
Crates per pallet	18
Qty per pallet	90
Cleaning	Use a mild soap solution to clean the seat. Please take care that the seat and hinges are not left damp. Dry off any residual water or cleaning agent with a soft cloth. Avoid contact of the seat and hinges with abrasive, corrosive or chlorine-based cleaners, as these may cause damage or lead to flash rust. Therefore, when cleaning the bowl, make sure that the seat and cover are kept in an upright position until all cleaning agent has been flushed away.



PRESSALIT®

Pressalit Style D300 1204

**Toilet seat with soft close and lift-off incl. hinge in stainless steel
Incl. design emblems, brushed brass and brushed steel**

Warranty

Your Pressalit toilet seat is subject to a 10-year product warranty against presentation of receipt. You should contact the dealer or store where you purchased the seat in the event of a complaint. In the event of complaints about defective materials and/or manufacturing defects and/or functional limitations ascertained within a 10-year period, following proper use and cleaning of the toilet seat (we reserve the right to examine the toilet seat at our factory), Pressalit will deliver spare parts or replace the toilet seat by an identical or similar model. If the color is no longer available, Pressalit reserves its right to replace the toilet seat by a similar model in the color white. Upon replacement or delivery of a spare part, the warranty period will not be extended nor renewed. If replacement is not possible, Pressalit may instead refund the purchase price.

Subject to changes in the product range, technical changes and printing errors.

pressalit.com

PRESSALIT®