




PRESSALIT®

Product Item number	Pressalit Public Classic Soap dispenser with sensor, 850 ml, brushed steel PY30715
Product description	VVS: 774565215 NRF: 1438639 RSK: 8920389 EAN: 5708590392738  Soap dispenser with stainless steel sensor for wall mounting. Drip-free soap pump with 850 ml container. Designed for regular liquid soap or hand cream. Lockable solution. Incl. mounting kit.
Materials	European stainless steel (AISI 316L)
Warranty Cleaning	A 2-year warranty is provided for material and manufacturing defects. The product is to be cleaned with hot water or a mild soapy cleaning agent. Do not use chlorine based or corrosive cleaners as it may result in the surface being damaged. Subject to changes in the product range, technical changes and printing errors.

pressalit.com