




# PRESSALIT®

<b>Product</b> <b>Item number</b>	Pressalit Public Classic Disinfection dispenser with sensor, 950 ml, incl. stand, brushed steel  PY30735
<b>Product description</b>	<p>VVS: 774565815 NRF: 1438618 RSK: 8920392 EAN: 5708590392349</p>  <p>Disinfection dispenser in brushed stainless steel. Mounted on floor standing stand. Disinfection pump with sensor and 950 ml tank. Designed for standard hand sanitizers. Lockable solution. Incl. mounting kit.</p> <p>Inspection opening indicates filling level.</p>
<b>Materials</b>	European stainless steel (AISI 316L)
<b>Warranty</b> <b>Cleaning</b>	<p>A 2-year warranty is provided for material and manufacturing defects.</p> <p>The product is to be cleaned with hot water or a mild soapy cleaning agent. Do not use chlorine based or corrosive cleaners as it may result in the surface being damaged.</p> <p>Subject to changes in the product range, technical changes and printing errors.</p>

[pressalit.com](http://pressalit.com)